

STAINLESS TURNTABLE SHEET “TS-SUS1” USER’S MANUAL


Thank you for purchasing the ViV laboratory’s product.

Please read this manual before using it.

1. CAUTION

- Be careful not to injure or cut your finger by the edge of “TS-SUS1”.
- The weight of the TS-SUS1 is approximately 1.7kg.
Please test and confirm that there is no problem with your turntable.

2. Setting



- Remove the old turntable sheet.
- Set the TS-SUS1 sheet as three-ball-point side coming down.
- Set a record disk on the TS-SUS1 sheet and put the attached disk stabilizer on it.
You can use your favorite disk stabilizer instead of attached one.
- Set the height of the tonearm if needed.

3. Specification

- dimension : ϕ 300 x 9 mm (max)
- weight : 1.7 kg (sheet)
- material : SUS304 stainless steel

ViV laboratory Ltd.

5-10-16 Imaizumidai, Kamakura city, Kanagawa pref., JAPAN 247-0045

www.vivaudiolab.com info@vivaudiolab.com